






## AXO DRESDEN GmbH

### - Applied X-ray Optics and High Precision Deposition -


- Development and production of multilayer X-ray optics
- Complex 1- and 2-dimensional X-ray optical systems
- Synchrotron optics from EUV to hard X-rays
- Monochromators for XRF
- High precision deposition by means of magnetron sputter deposition (MSD), large area pulsed laser deposition (LA-PLD) and dual ion beam deposition (DIBD) technologies
- Application and teaching in X-ray diffraction (XRD), X-ray reflectometry (XRR) and X-ray fluorescence (XRF)





#### HIGH PRECISION DEPOSITION



#### X-RAY OPTICAL SYSTEMS



#### MULTILAYER X-RAY OPTICS




# High Precision Deposition Techniques:

## Magnetron Sputter Deposition (MSD)

## Large Area Pulsed Laser Deposition (LA-PLD)


## Dual Ion Beam Deposition (DIBD)




Dual beam large area pulsed laser deposition plasma plumes caused by the interaction of focussed ns-laser pulses and the target surface


- High precision deposition of nanometer single and multilayers
- Accuracy and homogeneity in the picometer range
- Multilayer stacks for a wide variety of applications
- Constant and graded thickness distributions
- Deposition on rectangular or round substrates (typical dimensions: 8 inch diameter or 500 mm length)

## Application of nanometer single und multilayers for X-ray optics





## Parallel Beam X-ray Optics



Generation of a monochromatic parallel beam in one dimension

Parallel beam X-ray optics are optical components with a graded multilayer deposited on a substrate having a parabolic shape in beam direction. These optical devices convert (in one dimension) a divergent incoming beam into a parallel one, or vice-versa an incoming parallel beam into a focusing one.

In order to obtain high efficiency, the d-spacing of the multilayer has to be varied from the front end to the rear end of the optics in correspondence to the aspheric curvature. Either the X-ray source or the detector (or detector slit) may be placed at the optics focal distance, for primary or secondary side applications, respectively.



High precision 60 mm parallel beam optics on prefigured substrate (right) and on flat substrate, which are glued and bend after deposition (left)


Spectral lines	Cr, Co, Cu, Ga, Mo, Ag (others on request)
Mean Reflectivity	R > 70 %
Monochromacy	K $\alpha_1$ + K $\alpha_2$ or K $\beta$
Divergence	$\Delta\phi \leq 0.03^\circ$ (40 $\mu\text{m}$ source width)
Mirror length	L = 40...100 mm (on customer's request)
X-ray source geometry	line focus
Parallel beam width b	dependent on mirror length, geometry and X-ray wavelength
Typical b values	1.5 mm (Cu-K, L = 60 mm) 1.0 mm (Mo-K, L = 100 mm)
Geometry	typical focal length (source - mirror center) L = 60 mm, x <sub>m</sub> = 100 mm



## Convex Parallel Beam Optics




Convex curved parallel beam multilayer optics of various geometries for the generation of a compressed monochromatic parallel beam in one dimension



### Outstanding Feature

Increase of the parallel beam photon flux in combination with a focusing optic by reduction of beam width  $b$  compared to conventional parallel beam optics

Spectral lines	Cr, Co, Cu, Ga, Mo, Ag (others on request)
Mean reflectivity	$R > 70\%$
Monochromacy	$K\alpha_1 + K\alpha_2$ or $K\beta$
Mirror length	typical values: 20 mm – 40 mm (others on request)
X-ray source geometry	line focus
Geometry	customized



## Focusing X-ray Optics




Parallel beam and focusing X-ray optics with various geometries

**Focusing X-ray optics** are artificial optical components with a 1-dimensional lattice deposited on a substrate.

These optical devices convert a divergent incoming beam into a focusing one. To obtain high efficiency, the d-spacing of the lattice has to be changed from the front end to the rear end of the optics.

The device needs to have an elliptical figure of curvature to produce a focusing beam. The focus of the X-ray source is located in one of the two focal points of the ellipse.



Generation of a monochromatic focusing beam in one dimension

**Spectral lines** Cr, Co, Cu, Ga, Mo, Ag (others on request)

**Mean reflectivity** R > 70 %


**Monochromacy**  $K\alpha_1 + K\alpha_2$  or  $K\beta$

**Mirror length** typical values:  
L = 40 mm ... L = 80 mm  
(on customer's request)


**X-ray source geometry** line focus

**Focal line width b** dependent on spectral line, geometry and mirror length

**Geometry** customized




# Monochromators for XRF and Synchrotron Applications




Monochromatization of X-rays by means of a nanometer multilayer stack

Monochromators are optical devices with a 1-dimensional multilayer deposited on a substrate. To obtain high efficiency, the d-spacing of the multilayer is constant from the front end to the rear end of the monochromator. A plane figure is required for the monochromator. Depending on the application either high resolution or high flux multilayer monochromators can be fabricated.



Energy resolution and divergence of hard x-ray multilayer monochromators (\*)



Reflectivity of a Ni/C multilayer ( $d = 4.0 \text{ nm}$ ) in the spectral range between 8 keV ... 70 keV  
(Measurement: G. Falkenberg, Hasylab at DESY)


Spectral range	< 50 eV – 100 keV
Material systems	optimized on wavelength or on customer's request
Typical sizes	500 mm in length or 8 inches diameter
Resolution	$0.25 \% < \Delta E/E < 2 \%$ (periodic multilayers) $\Delta E/E > 5\%$ on request (aperiodic multilayers)
Thickness homogeneity	$\Delta d/d < 0.02 \%$
Applications	<b>monochromators</b> for laboratory X-ray sources and for synchrotrons, optimized for high reflectivity or tailored resolution <b>polarizers</b> in the soft X-ray range (O-K, Fe-L, Ni-L)

(\*) Ch. Morawe et al.- ESRF – Grenoble (F)




# Monochromators for X-ray Fluorescence Analysis Analysis of Light Elements (Be...Cl)

**High resolution multilayer  
Monochromators ( $2d = 3.02 \text{ nm}$ )  
vs. TIAP crystals ( $2d = 2.59 \text{ nm}$ )**



Performance of high resolution  
W/B<sub>4</sub>C multilayers and TIAP crystals  
for the detection of Si, Al and Mg

	<b>W/B4C</b>	<b>TIAP</b>
<b>Si K<math>\alpha</math>1,2</b>		
peak intensity [kcps]	<b>31.3</b>	<b>39.1</b>
P/B ratio	<b>55.6</b>	<b>129.2</b>
<b>Al K<math>\alpha</math>1,2</b>		
peak intensity [cps]	<b>37.6</b>	<b>36.8</b>
P/B ratio	<b>58.4</b>	<b>213.7</b>
<b>Mg K<math>\alpha</math>1,2</b>		
peak intensity [cps]	<b>107.8</b>	<b>103.5</b>
P/B ratio	<b>65.0</b>	<b>207.0</b>




Performance of high resolution  
multilayers (IWS30)\* at the Si K $\alpha$ -Line  
in comparison to W/B<sub>4</sub>C multilayers and  
TIAP (Thallium biphtalate) crystals.

	<b>W/B4C</b>	<b>IWS30</b>	<b>TIAP</b>
d [nm]	1.51	1.55	1.295
Order	1st	1st	1st
E [keV]	1.74	1.74	1.74
R [%]	21	26	ca. 25
E/ $\Delta E$	309	404	ca. 400


(\* ) in cooperation with Fraunhofer IWS, Dresden



## ASTIX-c 2-dimensional Parallel Beam X-ray Optics



ASTIX-c: collimating geometry




Parallel beam profile for  
Mo K $\alpha$  radiation  
 $I > 10^7$  cps  
(low power  $\mu$ -source)  
 $b^2 \approx 1 \text{ mm}^2$



**Collimating 2-dimensional X-ray optics  
in a modified Montel geometry<sup>(1)</sup>  
for the generation of 2-dimensional  
high intense parallel X-ray beams**

- Typical length 60 mm - 150 mm
- Application with all typical types of X-ray sources (rotating and fixed anodes, liquid metal jet and micro focus X-ray tubes)
- Typical parallel beam width:  
 $1 \text{ mm}^2 \leq b^2 \leq 5 \text{ mm}^2$
- Wavelengths: Cr, Co, Cu, Ga, Mo, Ag...
- High precision vacuum mirror housing




High precision vacuum mirror housings for  
ASTIX optics

(1) M. Montel - „The X-Ray Microscope with Catamegonic Roof-Shaped Objective“ in: *X-ray Microscopy and Microradiography*, Vol.5, 1957, pp 177 - 185



## ASTIX-f 2-dimensional Focusing X-ray Optics



ASTIX-f: focusing geometry

Focusing 2-dim. X-ray optics in a modified Montel geometry<sup>(1)</sup> for the generation of 2-dimensional high intense focused X-ray beams

- Typical length 60 mm - 150 mm
- Application with all typical types of X-ray sources (rotating and fixed anodes, liquid metal jet and micro focus X-ray tubes)
- Typical spot diameter <30 µm ... 500 µm
- Convergence: customized
- Wavelengths: Cr, Co, Cu, Ga, Mo, Ag, ...
- High precision vacuum mirror housing


### High Flux Optics

- high flux (HF) at sample position
- high integrated pixel intensity


### High Resolution Optics

- small spot size at sample position
- high resolution (HR) in the detector plane

Profile of diffracted beam between optics and focal point ( $f_2 = 310$  mm)



(1) M. Montel - „The X-Ray Microscope with Catameric Roof-Shaped Objective“ in: *X-ray Microscopy and Microradiography*, Vol.5, 1957, pp 177 - 185




# ASTIX-f

## Application of Different Types of Optics to a Lysozyme Sample

### High Resolution Optics

#### Advantages of HR optics

Small spot size at sample position  
High resolution (HR) in detector plane






#### Lysozyme sample

Measured with HR optics



### High Flux Optics

#### Advantages of HF optics

High integrated pixel intensity  
High flux (HF) at sample position




Nonius X-ray generator  
Cu K $\alpha$  radiation  
Marresearch dtb and  
mar345 imaging plate detector






## Twin Mirror Arrangement



**Twin Mirror Arrangement** consisting of primary mirror with housing (left) and secondary mirror with housing, beam tube and detector slit holder (right).




Setup of Twin Mirror Arrangement (TMA)

Upgrades available for all common X-ray instruments (for example Cr, Co, Cu, Ga, Mo, Ag K $\alpha$  or K $\beta$ , W L $\alpha$  or L $\beta$  radiation)  
Geometries on customers' request

A new quality in in-house X-ray reflectometry

### Special features

- Easy and fast sample alignment




X-ray reflectometry: Independence of peak intensities and angular positions at three different stage heights

- No influence of sample displacement errors up to 200 μm on peak position and intensity
- Sample alignment within **10 seconds**
- Dynamic range of more than **7 orders** of magnitude (cross intensity  $I > 3 \cdot 10^9$  cps)
- Low divergence ( Cu K $\alpha$  :  $\Delta\phi < 0.03^\circ$  )
- Detectable thin film thicknesses between 2 nm and 270 nm




## Twin Mirror Arrangement




**Twin Mirror Arrangement** consisting of primary mirror with housing (left) and secondary mirror with housing, beam tube and detector slit holder (right).

### A new quality in X-ray diffractometry – secondary parallel beam optics



Better resolution than soller slits because of a more than two times lower angular acceptance ( $Cu K\alpha : \Delta\phi < 0.03^\circ$ )




Powder diffraction measurement in the angular range of the quartz triplet in parallel beam geometry (TMA) (a) and BRAGG-BRENTANO geometry (b)

- Increased S/N-ratio due to sample fluorescence suppression
- Transmission higher than 60%
- Fits best to primary parallel beam optics
- Parallel beam geometry:
  - simplified sample preparation
  - increased accuracy


Superior  $K\beta$  suppression:  
 $I(Cu K\alpha_1) : I(Cu K\beta) > 1.000.000 : 1$   
 $(\theta - 2\theta - \text{scan at Si (400) wafer})$



## Thin Film XRF Reference Samples




**X-ray fluorescence reference sample** with up to 7 different elements on a silicon nitride membrane (above), PEEK sample holder (below).



Mass deposition of the elements on the reference sample with corresponding fluorescence emission line energies.


Element	Mass (ng/mm <sup>2</sup> )	Emission Lines (eV)	
		K $\alpha$	L $\alpha$
Pb	7.61±0.96	85335	10541
La	11.01±0.62	33298	4649
Pd	1.8±1.0	21123	2838
Mo	1.32±0.40	17444	2293
Cu	2.84±0.35	8040	930
Fe	5.04±0.87	6401	747
Ca	19.31±1.10	3691	341
Si	Substrate	1740	




Energy spectra of three 7-element reference samples. The energy range from ~2 keV to ~40 keV is covered with peaks of comparable intensity.

### Advantages of Thin Film XRF reference samples:

- Absorption free standard: no matrix correction necessary
- Substrate thickness of 100 nm / 200 nm permits transmission measurements and leads to low background from the substrate
- Mass depositions in the range of ng/mm<sup>2</sup> (1-3 atomic layers) permit quantification without the need to interpolate from higher values
- Uncertainty ≤ 1 ng/mm<sup>2</sup> (1 atomic layer)
- Wide selection of non-overlapping XRF lines, exact calibration curve with many points over a large energy range
- Signal strength easily adjustable by thickness, similar intensity for all elements
- High degree of uniformity & homogeneity (better than 1% for the full sample area)
- Application for adjustment of confocal μ-XRF possible
- Wide range of available elements (standard and tailored compilations)



Lateral homogeneity measured by XRF μ-beam mapping (1.2x1.2 mm<sup>2</sup> with 2.8x12 µm<sup>2</sup> beam) of La and Cu, and XRF large area mapping (15x15 mm<sup>2</sup> with 0.8x0.4 mm<sup>2</sup> beam) of La and Cu (from left to right)



## Teaching Tool for X-ray Reflectometry

X-ray reflectometry (XRR) is an essential technique for non-destructive characterization of surfaces, interfaces and thin films in the nanometer range. The Dresden University of Technology (TUD) Institute of Structural Physics (Nanostructural Physics - Prof. Dirk C. Meyer) offers an advanced practical course in X-ray reflectometry.




The samples characterized in these experiments are now available as a complete teaching tool in cooperation of TUD with AXO DRESDEN GmbH and FhG IWS. The scope of delivery also comprises an elaborated introduction in the basic principles of X-ray reflectometry and a commented trail guidance together with some protocols of students of these practical course of Dresden University of Technology as examples.

The set consists of three silicon substrates of  $40 \times 20 \text{ mm}^2$  size, one without coating, one with a thin single Ni-layer ( $d \approx 30 \text{ nm}$ ) and one with a 10-period Ni/C-multilayer (mean period thickness  $d \approx 4 \text{ nm}$ ). The films are deposited by Pulsed Laser Deposition (PLD).

The specular reflected intensity is measured as a function of grazing incident angle. Each measurement shows a characteristic pattern. Basic principles of total external reflexion and the interference of reflected and diffracted X-rays in a layer stack can be discussed. Layer thicknesses and interface roughnesses are determined by selected simulation programs that are on a CD delivered with the sample set.



Concept of sample selection:  
one interface bulk material / air  
↓  
single layer with two interfaces  
↓  
multilayer system with many interfaces





Simulation programs for quantitative analysis of experimental data are available as freeware [1]. Simulations of presented measurement were carried out using the program *IMD*.

[1] D.L. Windt, *IMD - Software for modeling the optical properties of multilayer films*, Computers in physics 12, 360-370 (1998).

### More information:

<http://www.physik.tu-dresden.de/isp/nano/lehre.php>



AXO DRESDEN GmbH - Applied X-ray Optics and High Precision Deposition - was founded as a Spin-off by employees of the Fraunhofer Institute Material and Beam Technology Dresden (IWS) with a participating share of the Fraunhofer Society.

Our production program contains both single X-ray optics and complex X-ray optical systems to generate high intensity collimated or focussed monochromatic X-ray beams, special customized depositions and applications in X-ray reflectometry and diffraction.


A wide assortment of flat and curved X-ray optics are available that can be used with the common types of X-ray sources with respect to wavelength, construction and focus geometry.



AXO DRESDEN develops and applies various high precision deposition techniques for the production of nanometer single and multilayers showing sub-nanometer precision.

To deliver our customers an optimum high quality solution we are working in close collaboration with the Fraunhofer IWS, further Fraunhofer Institutes, the University of Technology Dresden as well as with other national and international research institutions.

More than 18 years of experience in the fields of high precision deposition and design, development and application of multilayer X-ray optics guarantees our clientele effective and customized solutions in the fields of X-ray optics and high precision deposition.

We are looking forward to answer your questions and inquiries.





Die AXO DRESDEN GmbH ist ein Spin-Off aus dem Fraunhofer Institut für Werkstoff- und Strahltechnik (IWS) Dresden unter Beteiligung ehemaliger IWS-Mitarbeiter und der Fraunhofer Gesellschaft.

Unser Produktspektrum umfasst sowohl röntgenoptische Einzelkomponenten und komplexe strahlformende Systeme auf der Grundlage von Nanometer-Einzel- und Multischichten als auch spezielle Einzelbeschichtungen entsprechend Kundenwunsch sowie Applikationen im Bereich der Röntgen-Diffraktometrie und Reflektometrie.

Ein Standardsortiment an ebenen und unterschiedlich gekrümmten Röntgenoptiken für den Einsatz an den gebräuchlichen Röntgenröhren hinsichtlich Wellenlänge, Bauart und Fokusgeometrie steht zur Verfügung.


AXO DRESDEN entwickelt und benutzt unterschiedliche, sich ergänzende Präzisionsbeschichtungsverfahren, mit

denen Einzel- und Multischichten im Nanometerbereich mit sub-Nanometer-Präzision hergestellt werden können.

Um unseren Kunden stets optimale Lösungen anbieten zu können, unterhalten wir enge Kooperationsbeziehungen zum Fraunhofer IWS, zu anderen Fraunhofer Instituten, zur Technischen Universität Dresden sowie zu weiteren nationalen und internationalen Forschungseinrichtungen.

Unsere mehr als 18-jährige Erfahrung sowohl in der Präzisionsbeschichtung als auch auf röntgenoptischem Gebiet gewährleistet einem breiten Kundenkreis effektive, maßgeschneiderte Lösungen.

Wir freuen uns auf Ihre Anfragen.













# Nanometer thin films as XRF reference samples

M. Krämer<sup>1</sup>, R. Dietsch<sup>1</sup>, D. Weißbach<sup>1</sup>, G. Falkenberg<sup>2</sup>, R. Simon<sup>3</sup>, U. Fittschen<sup>4,5</sup>, T. Krugmann<sup>5</sup>

<sup>1</sup> AXO DRESDEN GmbH, Siegfried-Rädle-Str. 31, 01809 Heidenau, Germany; contact@axo-dresden.de

<sup>2</sup> HASYLAB at DESY, Notkestr. 85, 22603 Hamburg, Germany

<sup>3</sup> Institute for Synchrotron Radiation, FZ Karlsruhe, Hermann-von-Helmholtz-Platz 1, 76344 Eggenstein, Germany

<sup>4</sup> Chemistry Division, Los Alamos National Laboratory, P.O. Box 1663 Mail Stop K484, Los Alamos, NM 87545, USA

<sup>5</sup> Institute for Applied Chemistry, University of Hamburg, Martin-Luther-King-Platz 6, 20146 Hamburg, Germany

## Introduction



### Advantages:

- Absorption free standard – no matrix correction needed
- Wide range of suitable elements
- Selection of non-overlapping X-ray fluorescence lines
- Signal strength easily adjustable by layer thickness, comparable intensity for all elements
- High degree of uniformity and homogeneity possible (PVD)


### Applications:

- Quality assurance in XRF set-ups: alignment, optimization, calibration, minimum detection limits
- Confocal setups: depth resolution, detector capillary calibration

## Fabrication


Multielement samples were fabricated by PVD methods such as Magnetron Sputtering (MSD), Ion Beam Sputtering (IBS) and Pulsed Laser Deposition (PLD).

Polymer foils and silicon nitride membranes were used as ultrathin substrates.




## Energy spectrum


Sum spectra of the reference samples were measured at various X-ray sources. It can be seen that comparable characteristic emission lines of six to seven elements are well separated in the energy scale.



Sum spectrum of sample S10 measured at ANKA FLUO beamline (25 keV, 0.1x0.1 mm<sup>2</sup>)



Sum spectrum of sample RF4 measured at an Eagle laboratory μ-XRF with Rh tube (20 kV) at Los Alamos National Laboratory



## Homogeneity


Radiation source: ANKA FLUO-Beamline, KIT, Karlsruhe, Germany (shown here)  
Beamline L, HASYLAB at DESY, Hamburg, Germany  
μ-XRF Eagle III, Los Alamos National Laboratory, USA

Radiation energy: 9.5 keV, 25 keV, 26 keV, 28 keV synchrotron; 20 kV laboratory

Scanning area: 0.3x0.3 mm<sup>2</sup> up to 15x15 mm<sup>2</sup>  
Probe/Beam size: μm<sup>2</sup> range up to mm<sup>2</sup> range

Elements scanned: Si, Ca, La, Fe, Cu, Pb, Mo, Pd (La and Cu shown here)

### Homogeneity: Deviation < 1% over the entire scanned area




## Concentrations

### Chemical Analysis (University of Hamburg)

Foils were digested by HNO<sub>3</sub> / H<sub>2</sub>O<sub>2</sub> and a Mars Microwave digestion device. Elements were determined using a Ciroc CCD ICP-OES instrument (Spectro).

### XRF Analysis (ANKA, Karlsruhe and HASYLAB/DESY, Hamburg)

Sum spectra were calculated for numerous area mappings. Elemental concentrations were determined by ab initio calculations with instrumental parameters.



Comparison of calculated deposition amounts (Aim) and measurement results of different analysis methods (XRF and ICP-OES) for three identical reference samples (SF1, SF2, S100). Most target values were achieved and reproducible.

## Results

- Successful production of multi-element reference samples.
- High degree of homogeneity and reproducibility.
- Low X-ray absorption (below 2% in standard SF1 for E > 3 keV).
- Large spectral range without line overlap, emission lines have similar intensities.
- Membrane substrates withstand high radiation dose.